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1. Tonhazy NE, White NG, Umbreit WW. A rapid method for the estimation of the glutamic-aspartic transaminase in tissues and its application to radiation sickness. *Arch Biochem* 1950 28: 36-42.
2. Reitman S, Frankel S. A colorimetric method for the determination of serum glutamic oxalacetic and glutamic pyruvic transaminases. *Am J Clin Pathol* 1957 28: 56-63.
3. Murray, RL. Alanine aminotransferase. *Clinical chemistry theory, analysis and correlation*, 2ª ed. Kaplan LA, Pesce AP, comps. St. Louis: The C.V. Mosby Company. 1989: 895-898.
4. Wróblewski F, LaDue JS. Serum glutamic-pyruvic transaminase in cardiac and hepatic disease. *Proc Soc Exp Biol Med* 1956; 91: 569-571.
5. Bergmeyer, HU, M Horder. IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase. *J Clin Chem Clin Biochem* 1980;18:521-534.
6. Howe PE. The use of sodium sulfate as the globulin precipitant in the determination of proteins in blood. *J Biol Chem* 1921;49:93-107.
7. Howe PE. The determination of proteins in blood — a micro method. *J Biol Chem* 1921;49:109-113.
8. Wolfson WQ, y otros. A rapid procedure for the estimation of total protein, true albumin, total globulin, alpha globulin, beta globulin and gamma globulin in 10 ml of serum. *Am J Clin Pathol* 1948;18:723-30.
9. Saifer A, Gerstenfeld S, Vacsler F. Photometric microdetermination of total serum globulins by means of a tryptophan reaction. *Clin Chem* 1961;7:626-36.
10. Saifer A, Marven T. The photometric microdetermination of serum total globulins with a tryptophan reaction: a modified procedure. *Clin Chem* 1966;12:414-17.
11. Gendler, SM. Albumin. *Clinical chemistry theory, analysis and correlation*, 2ª ed. Kaplan LA, Pesce AP, comps. St. Louis: The C.V. Mosby Company. 1989: 1029-1033.
12. Webster D, Bignell AHC, Attwood EC. An assessment of the suitability of bromocresol green for the determination of serum albumin. *Clin Chim Acta* 1974 53: 101-108.
13. Louderback A, MealeyEH, Taylor NA. A new dye-binding technic using bromocresol purple for determination of albumin in serum. *Clin Chem* 1968 14: 793-794. (Abstract)
14. Pinnell AE, Northam BE. New automated dye-binding method for serum albumin determination with bromocresol purple. *Clin Chem* 1978;24:80-86.
15. King, EJ and Armstrong AR. A convenient method for determining serum and bile phosphatase activity. *Can Med Assoc J* 1934 31: 376-381.
16. Kind PRN, King EJ. Estimation of plasma phosphatase by determination of hydrolysed phenol with amino-antipyrine. *J Clin Pathol* 1954 7: 322-326.
17. Ohmori, Y. Uber die Phosphomonoesterase. *Enzymologia* 1937; 4: 217-231.
18. Fujita, H. Uber die Mikrobestimmung der Blutphosphatase. *J Biochem, Japan* 1939 30: 69-87.

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19. Petittlerc C, Delisle M, Martel M. Mechanism of action of Mg²⁺ and Zn²⁺ on rat placental alkaline phosphatase. I. Studies on the soluble Zn²⁺ and Mg²⁺ alkaline phosphatase. *Can J Biochem* 1975 53: 1089-1100.
20. Tietz NW, y otros. A reference method for measurement of alkaline phosphatase activity in human serum. *Clin Chem* 1983;29:751-761
21. Bowers GN Jr, y otros. IFCC methods for the measurement of catalytic concentration of enzymes. Part I. General considerations concerning the determination of the catalytic concentration of an enzyme in the blood serum or plasma of man. *Clin Chim Acta* 1979 98: 163F-174F.
22. Karmen, A. A note on the spectrophotometric assay of glutamic-oxalacetic transaminase in human blood serum. *J Clin Invest* 1955 34: 131-133.
23. Bergmeyer HU, y otros. Provisional recommendations on IFCC methods for the measurement of catalytic concentrations of enzymes. Part 2. IFCC method for aspartate aminotransferase. *Clin Chem* 1977 23: 887-899.
24. Bergmeyer HU, y otros. Provisional recommendations on IFCC methods for the measurement of catalytic concentrations of enzymes. Part 2. IFCC method for aspartate aminotransferase. *Clin Chem* 1978;24:720-721
25. Kramer B y otros. A simple technique for the determination of calcium and magnesium in small amounts of serum. *J Biol Chem* 1921;47:475-481.
26. Clark EP y otros. A study of the Tisdall method for the determination of blood serum calcium with suggested modification. *J Biol Chem* 1925;63:461-464.
27. Katzman E y otros. The determination of serum calcium by titration with ceric sulfate. *J. Biol Chem* 1937;118:539-544.
28. Cali y otros. A reference method for the determination of total calcium in serum. *En: Selected Methods of Clinical Chemistry*, Vol 8. Cooper GR, comp. Washington, DC: American Association for Clinical Chemistry. 1977: 3-8.
29. Kessler G y otros. An automated procedure for the simultaneous determination of calcium and phosphorus. *Clin Chem* 1964;10:686-703.
30. Michaylova V y otros. Photometric determination of micro amounts of calcium with arsenazo III. *Anal Chim Acta* 1971; 53: 194-198.
31. Scarpa A y otros. Metallochromic indicators of ionized calcium *Ann NY Acad Sci* 1978;307:86-112.
32. Ono T y otros. A new enzymatic assay of chloride in serum. *Clin Chem* 1988;34:552-553
33. Knoll VE y otros. Spezifische Kreatininbetimmung Im Serum. *Z Klin Chemi Clin Biochem.* 1970;8:582-587.
34. Haeckel R y otros. Simplified determinations of the "true" creatinine concentration in serum and urine. *J Cklin Chem Clin Biochem.* 1980;18:385-394.
35. Moss GA y otros. Kinetic enzymatic method for determining serum creatinine. 1975;21:1422-1426.
36. Jaynes PK y otros. An Enzymatic, Rate-Rate Assay For Serum Creatinine With a Centrifugal Analyzer. 1982; 28:114-117.
37. Fossati P y otros. Enzymatic Creatinine Assay: A New Colorimetric Method Based on Hydrogen Peroxide Measurement. 1983;29:1494-1496.
38. Whelton A y otros. Nitrogen metabolites and renal function. *En: Tietz Textbook of Clinical Chemistry*, 2^a ed. Burtis CA, Ashwood ER, comps. Philadelphia: W.B. Saunders Company. 1994;1513-1575.
39. Folin O y otros. A system of blood analysis. *J Biol Chem.* 1919;38:81-110.
40. Somogyi M. A reagent for the copper-iodometric determination of very small amounts of sugar. *J Biol Chem.* 1937;117:771-776.
41. Nelson N y otros. A photometric adaptation of the Somogyi method for the determination of glucose. *J Biol.* 1944;153:375-380.
42. Kaplan LA. Glucose. *Clinical chemistry theory, analysis and correlation*, 2^a ed. Kaplan LA, AJ Pesce AJ, comps. St. Louis: The C.V. Mosby Company. 1989: 850-856.
43. Berry MN y otros. Enzymatic determination of potassium in serum. *Clin Chem* 1989;35:817-20.
44. Van Pelt J. Enzymatic determination of sodium, potassium and chloride in serum compared with determination by flame photometry, coulometry and ion selective electrodes. *Clin Chem* 1994;40:846-847
45. Hubl W y otros. Enzymatic determination of sodium, potassium and chloride in abnormal (hemolyzed, icteric, lipemic, paraproteinemic, or uremic) serum samples compared with indirect determination with ion selective electrodes. *Clin Chem* 1994;40:1528-1531
46. Helgerson RC y otros. Host-guest Complexation. 50. Potassium and sodium selective-selective chromogenic ionophores. *J Amer Chem Soc* 1989;111:6339-50.
47. Kumar A y otros. Chromogenic ionophere-based methods for spectrophotometric assay of sodium and potassium in serum and plasma. *Clin Chem* 1988;34:1709-1712
48. Berry MN y otros. Enzymatic determination of sodium in serum. *Clin Chem* 1988;34:2295-2298
49. Malloy, HT y Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. *J Biol Chem* 1937; 119: 481-490.

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50. Meites, S. Bilirubin, direct reacting and total, modified Malloy-Evelyn method. En: Selected Methods of Clinical Chemistry, Vol. 9. Faulkner WR, Meites S, comps. Washington, DC: American Association for Clinical Chemistry. 1982: 119-124.
51. Murao S, Tanaka N. A new enzyme "bilirubin oxidase" produced by *Myrothecium verrucaria* MT-1. *Agric Biol Chem* 1981; 45: 2383-2384.
52. Osaki S, Anderson S. Enzymatic determination of bilirubin. *Clin Chem* 1984; 30: 971. (Abstract)
53. Perry B, y otros.. Measurement of total bilirubin by use of bilirubin oxidase. *Clin Chem* 1986; 32: 329-332.
54. Skeggs LT Jr. An automatic method for the determination of carbon dioxide in blood plasma. *Am J. Clin Pathol* 1960;33:181-185.
55. Korzun WJ, Miller WG. Carbon Dioxide. *Clinical chemistry theory, analysis and correlation*, 2^a ed. Kaplan LA, Pesce AP, comps. St. Louis: The CV Mosby Company. 1989: 869-872.
56. Koller A, Kaplan LA. Total serum protein. *Clinical chemistry theory, analysis and correlation*, 2^a ed. Kaplan LA, Pesce AP, comps. St. Louis: The C.V. Mosby Company. 1989: 1057-1060.
57. Reigler, E. Eine kolorimetrische Bestimmungsmethode des Eiweisses. *Z Anal Chem* 1914; 53: 242-245.
58. Weichselbaum, TE. An accurate and rapid method for the determination of proteins in small amounts of blood serum and plasma. *Am J Clin Path* 1946; 16: 40-49.
59. Dumas BT y otros. A candidate reference method for determination of total protein in serum. I. Development and validation. *Clin Chem* 1981; 27: 1642-1650.
60. Fales FW. Urea in serum, direct diacetyl monoxime method. *En: Selected Methods of Clinical Chemistry, Vol 9.* Faulkner WR, Meites S, comps. Washington, DC: American Association for Clinical Chemistry. 1982: 365-373.
61. Van Slyke y otros. A permanent preparation of urease, and its use in the determination of urea. *J Biol Chem*, 1914; 19: 11-228.
62. Fawcett JK y otros. A rapid and precise method for the determination of urea. *J Clin Pathol*, 1960;13:156-159.
63. Chaney y otros. Urea and ammonia determinations. *Clin Chem*, 1962;8:130-132.
64. Talke H y otros. Enzymatische Harnstoffbestimmung in Blut und Serum im optischen Test nach Warburg. *Klin Wochensh*, 1965;43:174-175.
65. Hallett y otros. Reduced nicotinamide adenine dinucleotide-coupled reaction for emergency blood urea estimation. *Clin Chim Acta*, 1971;35:33-37.
66. Patton y otros. Spectrophotometric and kinetics investigation of the Berthelot reaction for the determination of ammonia. *Anal Chem*, 1977;49:464-469.
67. Sampson EJ y otros. A coupled-enzyme equilibrium method for the measuring urea in serum: optimization and evaluation of the AACC study group on Urea Candidate reference method. *Clin Chem*, 1980;26:816-826.
68. National Committee For Clinical Laboratory Standards. Physician's office laboratory guidelines, tentative guideline, 2^a ed. NCCLS Document POL1-T2. Wayne, PA NCCLS, 1992.
69. Scott, M.G. *Electrolytes and Blood Gases. Tietz Textbook of Clinical Chemistry.* 3^a ed. Burtis CA, Ashwood ER, comps. Philadelphia: WB Saunders Company. 1999: 1058-1059.
70. National Committee for Clinical Laboratory Standards. Procedures for the handling and processing of blood specimens; tentative standard. NCCLS Document H18-A2. Wayne, PA NCCLS, 1999.
71. Overfield CV, y otros. Glycosis: a re-evaluation of the effect on blood glucose. *Clin Chim Acta* 1972 39: 35-40.
72. Rehak NN, Chiang BT. Storage of whole blood: effect of temperature on the measured concentration of analytes in serum. *Clin Chem*.1988; 34-2111-4.
73. Sherwin JE, Obernolte R. Bilirubin. En: *Clinical Chemistry: Theory, Analysis, and Correlation*, Kaplan LA, Pesce AJ, comps. St. Louis: The C.V. Mosby Company. 1989: 1009-1015.
74. Henry RJ y otros. *Clinical Chemistry: Principles and Technics*, 2^a ed. New York: Harper and Row, 1974: 417-421; 1058-1059.
75. Scott, M.G. *Electrolytes and Blood Gases. Tietz Textbook of Clinical Chemistry.* 3^a ed. Burtis CA, Ashwood ER, comps. Philadelphia: WB Saunders Company. 1999: 1065-1066.
76. National Committee for Clinical Laboratory Standards. Interference testing in clinical chemistry; proposed guideline. NCCLS Document EP7-A. Wayne, PA NCCLS, 2002.
77. Young DS. *Effects of drugs on clinical laboratory tests*, 3^a ed. Washington, DC: AACC Press, 1990.
78. National Committee for Clinical Laboratory Standards. How to define and determine reference intervals in the clinical laboratory, approved guidelines, 2^a ed NCCLS Document C28-A2. Wayne, PA NCCLS, 2000.
79. Painter PC, Cope JY, Smith JI. Reference Information for the clinical laboratory. En: Burtis CA, Ashwood ER, comps. *Tietz Textbook of Clinical Chemistry.* 3^a ed. Philadelphia: WB Saunders Company, 1999 págs. 676-684.
80. National Committee for Clinical Laboratory Standards. Evaluation of precision performance of clinical chemistry devices; approved guideline. NCCLS Document EP5-A. Wayne, PA: NCCLS, 1999.
81. National Committee for Clinical Laboratory Standards. Quality management for unit-use testing; proposed guideline. NCCLS Document EP18-A. Wayne, PA NCCLS, 2002.

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82. National Committee for Clinical Laboratory Standards. Method comparison and bias estimation using patient samples; approved guideline. NCCLS Document EP9-A2. Wayne, PA NCCLS, 2002.

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13. Bibliografia

1. Tonhazy NE, White NG, Umbreit WW. A rapid method for the estimation of the glutamic-aspartic transaminase in tissues and its application to radiation sickness. Arch Biochem 1950; 28: 36-42.
2. Reitman S, Frankel S. A colorimetric method for the determination of serum glutamic oxalacetic and glutamic pyruvic transaminases. Am J Clin Pathol 1957; 28: 56-63.
3. Murray, RL. Alanine aminotransferase. In: Clinical Chemistry: Theory, Analysis, and Correlation, 2nd ed. Kaplan LA, Pesce AP, eds. St. Louis: The C.V. Mosby Company. 1989: 895-898.
4. Wróblewski F, LaDue JS. Serum glutamic-pyruvic transaminase in cardiac and hepatic disease. Proc Soc Exp Biol Med 1956; 91: 569-571.
5. Bergmeyer HU, M Horder. IFCC methods for the measurement of catalytic concentration of enzymes. Part 3, IFCC method for alanine aminotransferase. J Clin Chem Clin Biochem 1980;18:521-534.

13. Bibliografia (segue)

6. Howe PE. The use of sodium sulfate as the globulin precipitant in the determination of proteins in blood. *J Biol Chem* 1921;49:93-107.
7. Howe PE. The determination of proteins in blood — a micro method. *J Biol Chem* 1921;49:109-113.
8. Wolfson WQ, et al. A rapid procedure for the estimation of total protein, true albumina, total globulin, alpha globulin, beta globulin and gamma globulin in 10 ml of serum. *Am J Clin Pathol* 1948;18:723-730.
9. Saifer A, Gerstenfeld S and Vacsler F. Photometric microdetermination of total serum globulins by Medias of a tryptophan reaction. *Clin Chem* 1961;7:626-636.
10. Saifer A, Marven T. The photometric microdetermination of serum total globulins with a tryptophan reaction: a modified procedure. *Clin Chem* 1966;12:414-417.
11. Gendler, SM. Albumin. In: *Clinical Chemistry: Theory, Analysis, and Correlation*, 2nd ed. In: Kaplan LA, Pesce AJ, eds. St. Louis: The C.V. Mosby Company. 1989: 1029-1033.
12. Webster D, Bignell AHC, Attwood EC. An assessment of the suitability of bromocresol green for the determination of serum albumin. *Clin Chim Acta* 1974; 53: 101-108.
13. Louderback A, MealeyEH, Taylor NA. A new dye-binding technic using bromocresol purple for determination of albumin in serum. *Clin Chem* 1968; 14: 793-794. (Abstract)
14. Pinnell AE, Northam BE. New automated dye-binding method for serum albumin determination with bromocresol purple. *Clin Chem* 1978;24:80-86.
15. King, EJ and AR Armstrong. A convenient method for determining serum and bile phosphatase activity. *Can Med Assoc J* 1934; 31: 376-381.
16. Kind PRN, King EJ. Estimation of plasma phosphatase by determination of hydrolysed phenol with amino-antipyrine. *J Clin Pathol* 1954; 7: 322-326.
17. Ohmori, Y. Uber die Phosphomonoesterase. *Enzymologia* 1937; 4: 217-231.
18. Fujita, H. Uber die Mikrobestimmung der Blutphosphatase. *J Biochem, Japan* 1939; 30: 69-87.
19. Petitclerc C, Delisle M, Martel M. Mechanism of action of Mg²⁺ and Zn²⁺ on rat placental alkaline phosphatase. I. Studies on the soluble Zn²⁺ and Mg²⁺ alkaline phosphatase. *Can J Biochem* 1975; 53: 1089-1100.
20. Tietz NW, et al. A reference method for measurement of alkaline phosphatase activity in human serum. *Clin Chem* 1983;29:751-761.
21. Bowers GN Jr, et al. IFCC methods for the measurement of catalytic concentration of enzymes. Part I. General considerations concerning the determination of the catalytic concentration of an enzyme in the blood serum or plasma of man. *Clin Chim Acta* 1979; 98: 163F-174F.
22. Karmen, A. A note on the spectrophotometric assay of glutamic-oxalacetic transaminase in human blood serum. *J Clin Invest* 1955; 34: 131-133.
23. Bergmeyer HU, et al. Provisional recommendations on IFCC methods for the measurement of catalytic concentrations of enzymes. Part 2. IFCC method for aspartato transaminasi. *Clin Chem* 1977; 23: 887-899.
24. Bergmeyer HU, et al. Provisional recommendations on IFCC methods for the measurement of catalytic concentrations of enzymes. Part 2. Revised IFCC method for aspartato transaminasi. *Clin Chem* 1978;24:720-721
25. Kramer B, et al. A simple technique for the determination of calcium and magnesium in small amounts of serum. *J Biol Chem* 1921;47:475-481.
26. Clark EP, et al. A study of the Tisdall method for the determination of blood serum calcium with suggested modification. *J Biol Chem* 1925;63:461-464.
27. Katzman E, et al. The determination of serum calcium by titration with ceric sulfate. *J. Biol Chem* 1937;118:539-544.
28. Cali et al. A reference method for the determination of total calcium in serum. *In: Selected Methods of Clinical Chemistry*, Vol 8. Cooper GR, ed. Washington, DC: American Association for Clinical Chemistry; 1977: 3-8.
29. Kessler G, et al. An automated procedure for the simultaneous determination of calcium and phosphorus. *Clin Chem* 1964;10:686-703.
30. Michaylova V, et al. Photometric determination of micro amounts of calcium with arsenazo III. *Anal Chim Acta* 1971; 53: 194-198.
31. Scarpa A, et al. Metallochromic indicators of ionized calcium *Ann NY Acad Sci* 1978;307:86-112.
32. Ono T, et al. A new enzymatic assay of chloride in serum. *Clin Chem* 1988;34:552-553.
33. Knoll VE, et al. Spezifische Kreatininbetimmung Im Serum. *Z Klin Chemi Clin Biochem.* 1970;8:582-587.
34. Haeckel R, et al. Simplified determinations of the “true” creatinine concentration in serum and urine. *J Cklin Chem Clin Biochem.* 1980;18:385-394.
35. Moss GA, et al. Kinetic Enzymatic Method For Determining Serum creatinine. 1975;21:1422-1426.
36. Jaynes PK, et al. An Enzymatic, Reaction-Rate Assay For Serum Creatinine With a Centrifugal Analyzer. 1982;28:114-117.
37. Fossati P, et al. Enzymatic Creatinine Assay: A New Colorimetric Method Based on Hydrogen Peroxide Measurement. 1983;29:1494-1496.

13. Bibliografia (segue)

38. Whelton A, et al. Nitrogen metabolites and renal function. In: Tietz Textbook of Clinical Chemistry, 2nd Ed. Burtis CA, Ashwood ER, eds. Philadelphia: W.B. Saunders Company. 1994;1513-1575.
39. Folin O, et al. A system of blood analysis. J Biol Chem. 1919;38:81-110.
40. Somogyi M. A reagent for the copper-iodometric determination of very small amounts of sugar. J Biol Chem. 1937;117:771-776.
41. Nelson N, et al. A photometric adaption of the Somogyi method for the determination of glucosio. J Biol. 1944;153:375-380.
42. Kaplan LA. Glucosio In: Clinical Chemistry: Theory, Analysis, and Correlation, 2nd ed. Kaplan LA, AJ Pesce AJ, eds. St. Louis: The C.V. Mosby Company. 1989: 850-856.
43. Berry MN, et al. Enzymatic determination of potassium in serum. Clin Chem 1989;35:817-820.
44. Van Pelt J. Enzymatic determination of sodium, potassium and chloride in serum compared with determination by flame photometry, coulometry and ion selective electrodes. Clin Chem 1994;40:846-847.
45. Hubl W, et al. Enzymatic determination of sodium, potassium and chloride in abnormal (hemolyzed, icteric, lipemic, paraproteinemic, or uremic) serum samples compared with indirect determination with ion selective electrodes. Clin Chem 1994; 40:1528-1531.
46. Helgerson RC, et al. Host-guest Complexation. 50. Potassium and sodium ion-selective chromogenic ionophores. J Amer Chem Soc 1989;111:6339-6350.
47. Kumar A, et al. Chromogenic ionophere-based methods for spectrophotometric assay of sodium and potassium in serum and plasma. Clin Chem 1988;34:1709-1712.
48. Berry MN, et al. Enzymatic determination of sodium in serum. Clin Chem 1988;34:2295-2298.
49. Malloy, HT and Evelyn KA . The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937; 119: 481-490.
50. Meites S. Bilirubin, direct reacting and total, modified Malloy-Evelyn method. In: Selected Methods of Clinical Chemistry, Vol 8. 9. Faulkner WR, Meites S, eds. Washington, DC: American Association for Clinical Chemistry; 1982: 119-124.
51. Murao S, Tanaka N. A new enzyme "bilirubin oxidase" produced by *Myrothecium verrucaria* MT-1. Agric Biol Chem 1981; 45: 2383-2384.
52. Osaki S, Anderson S. Enzymatic determination of bilirubin. Clin Chem 1984; 30: 971. (Abstract)
53. Perry B, et al. Measurement of total bilirubin by use of bilirubin oxidase. Clin Chem 1986; 32: 329-332.
54. Skeggs LT Jr. An automatic method for the determination of carbon dioxide in blood plasma. Am J. Clin Pathol 1960;33:181-185.
55. Korzun WJ, Miller WG. Carbon Dioxide. In: Clinical chemistry theory, analysis and correlation, 2nd ed. In: Kaplan LA, Pesce AJ, eds. St. Louis: The CV Mosby Company. 1989: 869-872.
56. Koller A, Kaplan LA. Total serum protein. In: Clinical Chemistry: Theory, Analysis, and Correlation, 2nd ed. In: Kaplan LA, Pesce AJ, eds. St. Louis: The C.V. Mosby Company. 1989: 1057-1060.
57. Reigler, E. Eine kolorimetrische Bestimmungsmethode des Eiweisses. Z Anal Chem 1914; 53: 242-245.
58. Weichselbaum, TE. An accurate and rapid method for the determination of proteins in small amounts of blood serum and plasma. Am J Clin Path 1946; 16: 40-49.
59. Dumas BT, et al. A candidate reference method for determination of total protein in serum. I. Development and validation. Clin Chem 1981; 27: 1642-1650.
60. Fales FW. Urea in serum, direct diacetyl monoxime method. In: Selected Methods of Clinical Chemistry, Vol 9. Faulkner WR, Meites S, eds. Washington, DC: American Association for Clinical Chemistry. 1982: 365-373.
61. Van Slyke, et al. A permanent preparation of urease, and its use in the determination of urea. J Biol Chem, 1914; 19: 11-228.
62. Fawcett JK, et al. A rapid and Precise method for the determination of urea. J Clin Pathol, 1960;13:156-159.
63. Chaney, et al. Urea and ammonia determinations. Clin Chem, 1962;8:130-132.
64. Talke H, et al. Enzymatische Harnstoffbestimmung in Blut und Serum im optischen Test nach Warburg. Klin Wochensch, 1965;43:174-175.
65. Hallett, et al. Reduced nicotinamide adenine dinucleotide-coupled reaction for emergency blood urea estimation. Clin Chim Acta, 1971;35:33-37.
66. Patton, et al. Spectrophotometric and kinetics investigation of the Berthelot reaction for the determination of ammonia. Anal Chem, 1977;49:464-469.
67. Sampson EJ, et al. A coupled-enzyme equilibrium method for the measuring urea in serum: optimization and evaluation of the AACC study group on Urea Candidate reference method. Clin Chem, 1980;26:816-826.
68. National Committee For Clinical Laboratory Standards. Physician's office laboratory guidelines, tentative guideline, 2nd ed. NCCLS Document POL1-T2. Wayne, PA: NCCLS, 1992.
69. Scott M.G. Electrolytes and Blood Gases. In: Tietz Textbook of Clinical Chemistry. 3rd ed. Burtis CA, Ashwood ER, eds. Philadelphia: WB Saunders Company. 1999: 1058-1059.

13. Bibliografia (segue)

70. National Committee For Clinical Laboratory Standards. Procedures for the handling and processing of blood specimens; tentative standard. NCCLS Document H18-A2. Wayne, PA: NCCLS, 1999.
71. Overfield CV, et al. Glycosis: a re-evaluation of the effect on blood glucose. *Clin Chim Acta* 1972; 39: 35-40.
72. Rehak NN, Chiang BT. Storage of whole blood: effect of temperature on the measured concentration of analytes in serum. *Clin Chem.*1988;34-2111-4.
73. Sherwin JE, Obernolte R. Bilirubin. In: *Clinical Chemistry: Theory, Analysis, and Correlation*, Kaplan LA, Pesce AJ, comps. St. Louis: The C.V. Mosby Company. 1989: 1009-1015.
74. Henry RJ, et al. *Clinical Chemistry: Principles and Technics*, 2° ed. New York: Harper and Row, 1974: 417-421; 1058-1059.
75. Scott M.G. Electrolytes and Blood Gases. In: *Tietz Textbook of Clinical Chemistry*. 3rd ed. Burtis CA, Ashwood ER, eds. Philadelphia: WB Saunders Company. 1999: 1065-1066.
76. National Committee For Clinical Laboratory Standards. Interference testing in clinical chemistry; proposed guideline. NCCLS Document EP7-A. Wayne, PA: NCCLS, 2002,
77. Young DS. *Effects of drugs on clinical laboratory tests*, 3° ed. Washington, DC: AACC Press, 1990.
78. National Committee For Clinical Laboratory Standards. How to define and determine reference intervals in the clinical laboratory, approved guidelines, 2nd ed. NCCLS Document C28-A2. Wayne, PA: NCCLS, 2000.
79. Painter PC, Cope JY, Smith JI. Reference Information for the clinical laboratory. In: Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: WB Saunders Company, 1999 pp. 676-684.
80. National Committee For Clinical Laboratory Standards. Evaluation of precision performance of clinical chemistry devices; approved guideline. NCCLS Document EP5-A. Wayne, PA: NCCLS, 1999.
81. National Committee For Clinical Laboratory Standards. Quality management for unit-use testing; proposed guideline. NCCLS Document EP18-A. Wayne, PA: NCCLS, 2002,
82. National Committee For Clinical Laboratory Standards. Method comparison and bias estimation using patient samples; approved guideline. NCCLS Document EP9-A2. Wayne, PA: NCCLS, 2002.

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© 2003, Abaxis, Inc., Unión City, CA 94587