



August 2008

VetCom

An Abaxis Newsletter For The Veterinary Community

Highlights

Special Look:

Potter Park Zoo

Amur Tiger Research Project
Utilizing the VetScan System

VetScan System in Galapagos

Monitoring Survival in Galapagos
Marine Iguanas

Case Study:

Critical Care Rotor

An Important Tool in the
Veterinary Practice

Spotlight:

August 20th Abaxis University
CE Course with Dr. Kent Adams

Authorized Distributors

Conference Calendar

Contact Us

Don't Miss A Single
Issue Of VetCom!



Welcome

Welcome to VetCom, our Free Newsletter available in print and online; providing valuable content to help you make the most of your VetScan® System.

Over 7,000 practices in North America similar to yours have invested in Abaxis VetScan products.

Thousands of practitioners are taking advantage of VetScan products because of its advantages over competitive in-house instruments including 20% lower operating costs, a 50% smaller sample size, 40% faster time to result, 10% less of hands-on time, sample integrity alerts, and on-board iQC® all assuring efficiency, quality and dependable results every time. Our commitment to you is non-stop development of new assays that address your needs at your clinic while NOW also offering 24/7 around the clock technical service.

The VetScan product line remains the most trusted in-house analyzers world-wide.

Did you know that our in-house laboratory analyzers are comparable to reference laboratory analyzers, and the sister product to the VetScan, the piccolo xpress® is FDA approved and used by Family Practitioners, Urgent Care Centers, Hospital Emergency Rooms, ICUs and CCUs. ■

Valerie Goodwin
Director, Marketing
Veterinary Medicine
valeriegoodwin@abaxis.com
(510) 675-6604



About Abaxis, Inc.

Abaxis is a leader in the development and commercialization of leading-edge innovative technology, tools and services that support best medical practices, enabling practitioners to respond to the health needs of their patients, while operating economical and profitable practices, at the point of care.

Abaxis develops, manufactures and markets portable blood analysis systems for use in any veterinary or human patient care setting to provide clinicians with rapid blood constituent measurements. The system consists of a compact, 5.1 kilogram (11.2 lbs), portable analyzer and a series of single-use plastic discs, called reagent discs that contain all the chemicals required to perform a panel of up to 13 tests on veterinary patients and 14 tests on human patients. The system can be operated with minimal training and performs multiple routine tests on whole blood, serum or plasma samples. The system provides test results in less than 12 minutes with the precision and accuracy equivalent to a clinical laboratory analyzer.

The veterinary business system, which consists of the VetScan VS2, the next generation chemistry, electrolyte, immunoassay, blood gas analyzer and the VetScan HM5 Hematology instrument, can be operated with minimal training and performs multiple routine tests on whole blood, serum or plasma samples.

Abaxis, Inc.
3240 Whipple Road
Union City, CA 94587
Tel. 800-822-2947
Fax 510-441-6150
www.abaxis.com

Subscribe today at
www.abaxis.com!



888-9300-REV P

August 2008

United States

American Veterinary Supply
800-869-2510

DVM Resources
877-828-1026

Great Western Animal Supply
800-888-7247

Hawaii Mega-Cor, Inc.
800-369-7711

IVESCO
800-457-0118

Lextron, Inc.
800-333-0853

Merritt Veterinary Supply
800-845-0411

Nelson Laboratories
800-843-3322

North East Veterinary Supply Co.
866-638-7265

Penn Vet Supply
800-233-0210

TW Medical
888-787-4483

Western Medical Supply
800-242-4415

Canada

Associated Vet Purchasing, Co.
604-856-2146

CDMV
450-771-2368

MidWest Drug
204-233-8155

VetNovations
705-726-9995

Vet Purchasing
519-284-1371

Vie et Sante
418-650-7888

Western Drug & Distribution
877-329-9332

Potter Park Zoo: Amur Tiger Research Project Using The VetScan System

Contributing Author: Tara Myers Harrison, DVM, MPVM, Veterinarian/Curator,
Potter Park Zoo, Lansing, MI 48912



Our project is quantifying the blood chemistry values of three, unanesthetized Amur tigers housed at our facility. This project is also part of a larger study involving operant conditioning (or training), and quantifying the normal levels of complete blood counts (CBC), cortisol levels, and blood pressures.

The eight subspecies of tiger (*Panthera tigris*) were once found throughout Asia. Today, only three subspecies remain within an ever-shrinking range of habitat. There are between 5,000 and 7,500 tigers left in the wild, but these populations are in danger from continued habitat loss and human interference. Maintenance of a stable, genetically viable captive population has become essential to this species' continued existence. The Species Survival Plan (SSP) in the United States is vital to this effort. Currently there are 151 Amur (*P. tigris altaica*), 60 Sumatran (*P. tigris sumatrae*), and 10 Indochinese (*P. tigris corbetti*) tigers under Tiger SSP management. The SSP makes breeding recommendations and produces a husbandry protocol. Knowledge of this species natural history and medical parameters is far from complete at this time.

Increasing numbers of institutions in the zoological setting are observing the benefits of conditioning and training to their animals, and incorporating it into their husbandry routines. There has been little evidence forth coming, however, that supports this on a physiological level. One recent study with chimpanzees (*Pan troglodytes*) found that positive reinforcement training significantly decreased the blood variables indicative of stress. The samples being tested were obtained without the use of anesthetics. The three tigers were born at Potter Park Zoo and have been trained to allow voluntary blood collection as well as several other veterinary procedures.

This study is gathering baseline data on the general chemistry, blood counts, blood pressure, and serum cortisol levels of unanesthetized tigers, animals that have been trained through positive reinforcement techniques. We aim to determine if differences in these variables can be observed when comparing anesthetized and unanesthetized animals and animals that are trained or are anesthetized without using hand-injection techniques. This data will assist veterinarians in better understanding how to evaluate blood work from an anesthetized large felid, improving the health of captive populations. ■

- Jul. 19-22 **AVMA**
New Orleans Hilton
New Orleans, LA
- Jul. 27-31 **World Veterinary Conf.**
Vancouver Convention
Center
www.vcea.ca
- Aug. 1 **Hambletonian**
Sheraton Meadowlands
East Rutherford, NJ
- Aug. 10 **VMSG**
Ventura, CA
www.vmsg.com
- Aug. 11-13 **Association of Avian
Veterinarians**
Hyatt Regency Hotel
Savannah, GA
www.aav.org
- Aug. 13-16 **South Dakota VMA**
Rapid City, SD
www.sdvmed.org
- Aug. 14-17 **Keystone Vet
Convention**
Hershey Lodge &
Convention Center
Hershey, PA
www.pavma.org
- Aug. 14-17 **Southern Veterinary
Conference**
Birmingham, AL
- Aug. 15-17 **XXI Annual Convention**
San Juan, Puerto Rico
www.cmvpr.org
- Aug. 23-26 **CVC Central**
Kansas City, MO
www.thecvc.com
- Sep. 17-19 **Purdue University**
School of Vet Med
Lafayette, IN
- Sep. 17-21 **IVECCS**
JW Marriott Desert
Ridge
Phoenix, AZ
www.veccs.org



This anesthetized Amur tiger is receiving a physical examination at Potter Park Zoo. Evaluation of blood results from this animal while awake may be different from when she is anesthetized.



This tiger is being trained through positive reinforcement while in a specialized cage to accept veterinary procedures such as drawing blood from the tail vein.

Potter Park Zoo Update

Potter Park Zoo in conjunction with Abaxis is evaluating the health of Amur tigers through operant conditioning or behavioral training. Blood is collected from the lateral tail vein of the tiger through the bars. Evaluating blood samples of animals trained for medical procedures is important for the evaluation of this endangered species. In the past, the only way to monitor the health of these animals is through anesthesia. Various anesthetic protocols and methods of inducing anesthesia in these animals could produce a stress response in various biochemical analytes as well as alter biochemical results. All of the normals developed for these species are based on animals that are anesthetized. Evaluating and developing normals for animals that are unanesthetized will allow veterinarians to determine the effect of anesthesia upon the biochemical analytes as well as to determine the health of an animal without the risk of anesthesia and also, in theory remove the stress response associated with anesthesia.

In addition to evaluating all the biochemical profiles on the Comprehensive Diagnostic Profile rotor, we are also evaluating cortisol levels to determine the stress of these animals that are trained for this procedure and exhibit no outward signs of stress.

Preliminarily, the ISIS (International Species Inventory System, Zoo Standard normals) mean normal for glucose is 134 mg/dL and our results are showing the overall mean for these animals is 100.4 mg/dL. Additional samples are being analyzed, as well as additional profiles, but at this point the trend is showing a much lower glucose level for awake, trained animals than animals that were anesthetized.

We will continue to analyze this and other analytes including cortisol, alkaline phosphatase, alanine amino transferase and creatinine. We hope that these results will aid in zoos and other conservation organizations to continue to provide improving health care for this endangered species. ■

Potter Park Zoo and its staff would like to thank Abaxis and Mr. Zovac for their participation and involvement in our research project involving Amur tigers.

1-800-822-2947

VetCom Issues

Valerie Goodwin
VetCom, EDITOR
Director, Veterinary Marketing
510-675-6604
valeriegoodwin@abaxis.com
questions and comments

Customer Service

David Bannister, ext. 6540
Director - Customer Care

Technical Support

Craig Darien, ext. 1521
Supervisor, Technical Support

North American Sales

Jeff Sumpter, ext. 1335
Area Sales Manager,
South Central
jeffsumpter@abaxis.com

Torrey Inez, ext. 1427
Area Sales Manager, South-
west
torreyinez@abaxis.com

Bill Wilson, ext. 1455
Area Sales Manager, Southeast
billwilson@abaxis.com

Andy Koupas, ext. 1424
Area Sales Manager, West
andykoupas@abaxis.com

John Therrien, ext 1476
Area Sales Manager,
North Central
johntherrien@abaxis.com

Gerard Cabrera, ext. 1448
Area Sales Manager, Northeast
gerardcabrera@abaxis.com

Chris Ballinger, ext. 1492
Area Sales Manager, Midwest
chrisballinger@abaxis.com

VetScan System In Galapagos: Monitoring Survival In Galapagos Marine Iguanas

Contributing Authors: L. Michael Romero (Tufts University), Martin Wikelski (Princeton University), Gregory L. Florant (Colorado State University)



Many wild animal populations experience periodic catastrophic population declines. There can be many different causes for these declines, such as new diseases (e.g. introduction of West Nile virus to North America), widespread famines (e.g. during droughts), and human-induced changes in environmental conditions (e.g. climate change or habitat alteration). Wildlife managers around the world are searching for factors that allow them to predict which individual animal will most likely be in distress and needs help during environmental changes.

Animals have a number of physiological systems that help them cope with these events, but it is very difficult to study which physiological responses are most important for survival and why those responses are insufficient for the majority of individuals. The main problem is that it is very hard to study wild animals that are under significant physiological stress and may be in the process of dying.

For the past few years we've been using the Galapagos Islands as a natural laboratory to study environmentally induced periods of stress that lead to mortality. The Galapagos are subject to periodic El Niño weather events. The main effect of El Niños on the marine environment of the Galapagos is the failure of nutrient-rich upwelling. The lack of nutrients turns one of the most productive marine environments in the world into a natural disaster. The bottom of the food web contracts, leaving all the higher trophic levels in deep trouble. Mortality strikes the fish, the sea birds, and the sea lions.

We have studied the Galapagos marine iguana (*Amblyrhynchus cristatus*) population on Santa Fe Island continuously since 1981 and witnessed massive starvation - up to 90% - during an El Niño which can last for several months in a bad year. The Galapagos marine iguanas are the only marine lizard. Apparently due to a reliance on a unique suite of endosymbiotic gut bacteria, these iguanas feed exclusively on intertidal and subtidal marine algae, which virtually disappears during an El Niño. We have a 24 year data set on body length and weight from known individuals and have shown that mortality during El Niños is different between different size and weight classes; the largest individuals suffer the highest mortality. Perhaps as a coping mechanism, many iguanas shrink in body



Medical & Research



Universities



Veterinary Facilities



length during El Niños. Starting in 1997, we have also been collecting physiological data on how individuals respond to long-term fasting which represents a unique stress in the tropics.

We have been making use of the fact that Galapagos marine iguanas repeatedly experience natural catastrophes (El Niños) but at the same time are tame and can easily be observed as they die. What we've been trying to understand is why some animals live and why others die. All the animals face identical environmental conditions, yet some animals will be resting along the coastline next to carcasses. Due to their tameness, marine iguanas uniquely allow us to repeatedly take blood samples and therefore we can determine the blood chemistry from individually-known animals with a known growth trajectory and life history throughout a population bottleneck.



The Vetscan has been an important part of our recent work. We are using the natural experiment of coping with El Niño conditions to understand whether different blood chemistry parameters correlate with survival. Our goal is to be able to predict before the El Niño which animals are likely to survive and which are likely to die. This is a critical question in ecology. It is currently unknown what characteristics allow an individual, particularly marine reptiles, to survive

environmental emergencies (e.g. long term fasting) while its neighbors succumb, even though this is an important foundation in natural selection.

Using the VetScan, we have already found some surprising results. For instance, all Galapagos marine iguanas have extremely high sodium levels. This is probably not surprising for a marine animal that drinks copious sea water as they forage on a marine algae. However, despite having a salt gland, the marine iguanas thrive with higher blood sodium concentrations than the Vetscan avian/reptilian rotors can measure. Much of the rest of our data are still being analyzed, but we are excited about the data the Vetscan has been able to provide.

Understanding how the blood chemistry changes in this reptile under various environmental and physiological conditions may very well enhance our ability to provide conservation measures to this and other marine animals. In addition, we hope that our data will provide us with clues about how natural selection works at the physiological level to promote changes in population structures. ■



Avian/Reptilian Profile Plus
12 tests per rotor

This rotor provides the following parameters: ALB, AST, BA, CA, CK, GLOB, GLU, K+, NA+, PHOS, TP, UA

An Important Tool in the Veterinary Practice

Contributing Author: Craig Tockman, DVM

The critical care rotor is often underutilized in the veterinary practices using the Vetscan point-of-care chemistry analyzer. In reality, it has incredible utility and you can and should make its utilization part of your everyday routine.

The testing parameters of the critical care rotor make it an excellent choice for monitoring your hospitalized patients on fluid therapy, regardless of the disease process. It is also an excellent rotor for recheck monitoring for many types of disease. A brief look at the tests the rotor offers will help understand how useful this rotor can be for your patients and your practice.

BUN and Creatinine allow the clinician to monitor renal function and hydration. ALT (alanine aminotransferase) maintains a window on liver cell death or leakage. Blood glucose measurements are reported. Sodium and potassium have always been present on the rotor, and now the addition of chloride provides the veterinarian with a complete electrolyte panel for fluid therapy determination.

We particularly like to monitor our chronic renal failure patients, Addisonian patients and diabetics with the critical care rotor. Simply monitoring BUN and Creatinine in renal patients is insufficient as these patients are often older or have other concurrent disease. The ability to monitor renal values along with glucose, ALT and electrolytes provides a better overall picture of the patients' health status between full diagnostic panels. Utilizing this rotor will help you catch early signs of hepatic disease (via ALT) and electrolytes are vital in renal patients, especially when the need for fluid therapy or diuretics is indicated. With routine use of this rotor, we have also found a significant number of renal patients with low levels of TCO₂ indicating a chronic metabolic acidosis. The patients tend to feel worse, have poor appetites and do not thrive on standard treatments. We are able to identify them early and correct the metabolic disorder more effectively.

Glucose monitoring alone in the diabetic is also insufficient for complete patient assessment. Diabetic cats often develop hepatic lipidosis and renal failure, making concurrent monitoring of BUN, Creatinine and ALT vital. The importance of electrolyte monitoring in these patients should not be underestimated. The need for monitoring renal and hepatic values, along with complete electrolytes is obvious for the Addisonian patient. In our hospital, the critical care rotor is used every 3-4

months to monitor these patients, along with a CBC. A comprehensive diagnostic rotor plus CBC is used every year.

Hospitalized patients on fluid therapy require frequent monitoring, regardless of the reason for the hospitalization. The critical care rotor provides an excellent panel to monitor these patients at least once daily. Possibly the most overlooked value is tCO₂ (total carbon dioxide). Complete blood gas analysis is not often practical in the veterinary practice, but there are times when it is needed.

The availability of tCO₂ on the critical care rotor provides you with a good idea as to the acid/base status of the patient. This, in combination with the electrolyte levels, allows you to choose the most appropriate and life saving fluid therapy.

When handled anaerobically, total carbon dioxide levels are equivalent to bicarbonate levels in the body. It is vital that blood drawn for this test be placed in the lithium heparin tube immediately after the sample is obtained to minimize contact with the air.

tCO₂ determination can help the clinician determine if more extensive blood gas analysis is needed while allowing refinement of the diagnostic and therapeutic plan. The combination of tCO₂, clinical signs and other blood values allows for a matrix that allows more accurate interpretation than one value alone. The following table summarizes the diagnostic value of tCO₂:

Total CO ₂ Level	Acid-Base Disorder
Decreased tCO ₂	Metabolic Acidosis Respiratory Alkalosis
Increased CO ₂	Metabolic Alkalosis Respiratory Acidosis

Since tCO₂ is equivalent to Bicarbonate levels, we can narrow the likely metabolic disorder based on its level.

First let's look at decreased Total Carbon Dioxide. Decreased tCO₂ levels indicate either metabolic acidosis (the most common acid-base disorder) or respiratory alkalosis. The most common causes of respiratory alkalosis are: tachypnea due to hypoxia, pulmonary disease, direct stimulation of the medullary respiratory center and excessive mechanical ventilation. Therefore, if the clinical signs of the patient are not consistent with these types of disorders, then metabolic acidosis is the likely reason for a decreased total carbon dioxide.

The most common causes of metabolic acidosis include:

- Small bowel diarrhea causing a loss of bicarbonate rich fluid. This can lead to hyperchloremic metabolic acidosis, so the chloride levels provided on the rotor can be very helpful in this determination.
- Ethylene glycol toxicity, salicylate toxicity and some medication can cause metabolic acidosis.
- Metabolic production of acid caused by lactic acidosis or diabetic ketoacidosis.
- Excessive urinary loss of bicarbonate due to renal failure, hypoadrenocorticism, etc.

A patient with a low tCO₂ can be assumed to be in a state of metabolic acidosis when seen with the disease states listed above.

What if the total carbon dioxide is elevated? Elevated total carbon dioxide indicates the possibility of metabolic alkalosis or respiratory acidosis. Metabolic alkalosis is most commonly caused by loss of chloride rich fluid from the intestinal tract or kidneys, and chronic alkali administration. Vomiting and chronic diuretic administration would be the most common causes of these conditions.

Respiratory acidosis is most commonly a result of hypoventilation due to airway obstruction, neuromuscular disease, neuromuscular blocking drugs, restrictive respiratory defects (diaphragmatic hernia, pleural effusion, etc.), pulmonary disease or inadequate mechanical ventilation.

Utilizing clinical signs or confirmed diagnostics together with the total carbon dioxide and chloride levels will help determine the most likely acid-base status. Based on this, you can determine the proper fluid therapy and drug administration for your patients.



Courtney, a 14 year-old spayed female Himalayan cat presented for PU/PD, excessive appetite and weight loss. 14 months prior to presentation, she was diagnosed with polycystic renal disease with a concurrent normochromic, normocytic anemia. The patient had been maintained on a renal diet and was receiving injection of Epoetin Alpha to maintain the red cell mass. Physical examination revealed a thin patient (BCS 4/9) and palpably enlarged and cystic kidneys. Initial CBC and chemistry results were as follows:

Case Study Utilizing the Critical Care Rotor

WBC: 7.25 (5.5-19.5)
Lym: 1.36 (1.5-7)
Ly%: 18.7 (20-55)
Mon: 0.39 (0-1.5)
Mo%: 5.4 (1-3)
Gra: 5.50 (2.5-14)
Gr%: 75.9 (35-80)
RBC: 5.24 (5-10)
HGB: 8.7 (8-15)
HCT: 23.21 (24-45)

MCV: 44 (39-55)
MCH: 16.5 (12.5-17.5)
MCHC: 37.3 (30-36)
RDWc: 19.2%
PLT: 396 (300-800)
BUN: 99 (10-30)
CRE: 4.1 (0.3-2.1)
ALP: 24 (10-90)
ALT: 48 (20-100)
TBIL: 0.3 (0.1-0.6)

AMY: 1630 (300-1100)
CA: 11.2 (8-11.8)
PHOS: 10.4 (3.4-8.5)
GLU: 86 (70-150)
NA+: 153 (142-164)
K+: 4.9 (3.7-5.8)
TP: 7.9 (5.4-8.2)
ALB: 3.8 (2.2-4.4)
GLOB: 4.0 (1.5-5.7)
T4: 1.8 (1.5-4.8)

The patient was placed on 0.9% NaCl at 1.5 times maintenance to begin diuresis for a 48 hour period. Two days after initial presentation, a CBC/critical care profile was performed. The anemia had worsened (RBC 4.37, HCT 19.70), but was still normochromic, normocytic. BUN and CRE were decreased (82 and 3.5), K+ was low normal at 3.7 (3.7-5.8) and total CO₂ was 14 (Normal 15-24) indicating a metabolic acidosis due to renal failure.

Fluids were changed to Lactated Ringers solution with the addition of sodium bicarbonate using the formula:

$$0.5 \times \text{Body Weight (in Kg)} \times (\text{Desired tCO}_2 \text{ mEq/L} - \text{Measured tCO}_2 \text{ mEq/L}) = \text{mEq Bicarbonate Required}$$

The following day, CBC and critical care panel results were:

WBC: 7.07 (5.5-19.5)
Lym: 0.83 (1.5-7)
Ly%: 11.7 (20-55)
Mon: 0.34 (0-1.5)
Mo%: 4.9 (1-3)
Gra: 5.89 (2.5-14)
Gr%: 83.4 (35-80)
RBC: 4.28 (5-10)
HGB: 6.7 (8-15)
HCT: 18.78 (24-45)
MCV: 44 (39-55)

MCH: 16.5 (12.5-17.5)
MCHC: 37.3 (30-36)
RDWc: 19.2%
PLT: 396 (300-800)
BUN: 73 (10-30)
CRE: 4.2 (0.3-2.1)
ALT: 66 (20-100)
GLU: 80 (70-150)
NA+: 141 (142-164)
K+: 3.8 (3.7-5.8)
tCO₂: 17 (15-24)



Valerie Goodwin
 Director Marketing - Veterinary
 Medicine
 VetCom Editor
 510-675-6604
 valeriegoodwin@abaxis.com



Craig M. Tockman, DVM
 Director of Professional
 Services, Abaxis, Inc.
 800-822-2947, ext. 1422
 craightockman@abaxis.com



Baerbel Koehler
 Abaxis Sales and Marketing
 Manager for Europe, Africa
 and the Middle East
 Darmstadt, Germany
 +49 6151 350 79 0
 baerbelkoehler@abaxis.de



Kent Adams, DVM
 Abaxis Advisory Board Member
 Equine/Large Animal Practitioner
 kentadams@abaxis.com



Don J. Harris, DVM
 Abaxis Advisory Board Member
 Avian/Exotic Practitioner
 donharris@abaxis.com

The fluid change and addition of bicarbonate had corrected the metabolic acidosis. The anemia had stabilized, but the sodium level was slightly low. Fluids were again changed to 0.9% NaCl at 1.5 times maintenance, with bicarbonate added and calculated to reach a tCO₂ of 19-20. After 2 days, the critical care rotor results were:

BUN: 68 (10-30) NA+: 146 (142-164)
 CRE: 4.2 (0.3-2.1) K+: 3.8 (3.7-5.8)
 ALT: 58 (20-100) tCO₂: 19 (15-24)
 GLU: 95 (70-150)

The patient was released to the owner with a guarded to poor prognosis, as the anemia continued to worsen. However, the cat remains stable with regular subcutaneous fluid therapy and Epopen three times weekly.

This case clearly shows the clinical utility of the critical care rotor. Correction of metabolic acidosis and monitoring electrolyte balance is an important part of fluid therapy and helps improve recovery. ■

Win A Free Upgrade!

Submit A Case Study And Win!

If you have a case study that includes how your in-clinic laboratory system has made a difference, and your case study is published in VetCom - you can WIN!



What you will WIN

You will be given a choice of a Free Upgrade!

You Choose:

- a. HMT or HM2 to the all new HM5
- b. VetScan Classic to VS2



How to enter

Send your case studies to Valerie Goodwin.

E-mail: valeriegoodwin@abaxis.com

Mail: Valerie Goodwin
 Marketing Director,
 Veterinary Medicine
 Abaxis, Inc.
 3240 Whipple Road
 Union City, CA 94587

(Open to owners and principals of veterinary practices, research facilities, academic environments, and pharmaceutical/biotech companies. Instrument being upgraded must be returned to Abaxis.)

VetScan HM5, a fully-automated 5-part cell counter specifically designed for veterinary applications, offers a comprehensive 22-parameter complete blood count with cellular histograms in just minutes. Its superior performance, elegant design, ease of use, true database management capability, and minimal maintenance make it the optimal hematology system for veterinary clinics, research laboratories, and pharmaceutical/biotech companies.

THREE EASY STEPS

Seconds of hands-on time are needed for performing CBC measurements using the VetScan HM5. All you need to do is:



1) enter information



2) insert sample and press start




3) read results

FEATURES

- 50 µl sample size
- 5-part differential on Cat, Dog, Horse
- Seven other species available in 3-part differential mode
- Small, space-saving, and elegant design with large LCD display
- Pre-dilute mode available for challenging samples
- Simple, intuitive, and easy to navigate software interface
- 22 parameters, including eosinophil counts and eosinophil percentage
- Advanced and integrated self-cleaning system
- Large database storage capacity, simplicity, flexibility, and usability
- CD drive and USB ports
- Integrates to the VetScan and data management systems
- User-reminder functions
- Multi-user mode

BENEFITS

- 10 species available now
- Small footprint
- Easy to use
- Reliable results, including out of range values
- Rapid test turn-around with minimal hands-on time
- Customizable reference ranges
- Optimized performance with minimal maintenance
- Automatically saves up to 2000 test records. All downloadable through USB pen drive and integrated data management system
- Allows for software updates and future instrument upgrades
- Integrated result reporting and comprehensive laboratory analysis
- Software tracks reagent usage and notifies user to replace reagents when necessary
- Eco-friendly reagents



COMBINED COMPREHENSIVE HEMATOLOGY AND CHEMISTRY REPORT

The test results and four histograms are displayed in minutes on an 8.5 x 11 printout.

Call 1-800-822-2947 or email vetscansales@abaxis.com to learn more about the HM5 (5-part differential) cell counter.

ABAXIS
UNIVERSITY



You are invited to a **COMPLIMENTARY** **1 HOUR ONLINE CE** **COURSE**


on August 20th, 2008
sponsored by **ABAXIS**

Designing and Implementing an Effective Equine Wellness Program

Program Outline:

- Demographic changes that make shift to equine wellness programs critical
- Using wellness programs to help establish brand identity
- Principles of a solid program
- Steps to successful implementation
- Expectations of marketing
- Reasons for success/failure
- Understanding the value of lab results
- Financial results of an aggressive program

Dr. Kent Adams
owner of Appalachian
Veterinary Services
in Riner, Virginia,
shares his approach
and success with
practitioners on how
to implement an
equine wellness
program.



**REGISTER NOW FOR AUGUST 20th at 11:30am OR 1:30pm (CST)
TO START YOUR COMPLIMENTARY
1 HOUR CE COURSE THROUGH
THE ABAXIS UNIVERSITY WEBSITE**



www.abaxisuniversity.com

or call us at 1.800.822.2947 ext. 6585

CE Courses meets the requirements for 1 hour of continuing education credit in jurisdictions which recognize AAVSB RACE approval; however participants should be aware that some boards have limitations on the number of hours accepted in certain categories and/or restrictions on certain methods of delivery of continuing education.

ABAXIS
innovation you can **trust**



Name: Harris Whitford
Doctor: Captain Robert Barhan, MD
Instrument: piccolo xpress

Name: Shadow
Doctor: Captain Grant Suski, DVM
Instrument: VetScan VS2

Millions of soldiers rely on Abaxis technology. **You can too.**

When staff sergeant Harris Whitford noticed an increasing pattern of sluggish behavior from his patrol dog, the two marched directly to the army veterinarian. It was no surprise when Shadow was diagnosed with hypothyroidism, but sergeant Whitford couldn't believe that little blue analyzer used to monitor Shadow's hormone levels was the same little blue analyzer the army used to test his blood during a routine physical. Strange coincidence? We don't think so. Abaxis instruments are rapidly becoming the standard of care for physicians and veterinarians around the world. We're committed to developing leading edge technology, tools and services that support best medical practices for all types of health care professionals. Our instruments are trusted by the Department of Defense, doctors' offices, emergency rooms, research labs and veterinary clinics everywhere for fast, reliable laboratory results. We believe that everyone is entitled to the highest quality health care available. And with our technology and your expertise we hope to make that possible for all the staff sergeants and Shadows of the world.

Call us today to schedule a live, on-site demonstration.

Tel: 510.675.6500 | 800.822.2947 | www.abaxis.com



Get the Most Out of Your VetScan Chemistry Analyzer

Comprehensive Diagnostic Profile

Ideal For: Providing a complete diagnostic analysis for pre-anesthetic testing, general health profile, ill patient diagnostics, geriatric testing, and wellness testing.

Analytes Include: ALB, ALP, ALT, AMY, BUN, CA, CRE, GLOB, GLU, K+, NA+, PHOS, TBIL, TP

Used For: Wellness testing, pre-anesthetic testing, critical cases, diagnostic testing and ill patient examinations, re-check examination for chronic disease.

Critical Care Plus

Ideal For: Serial testing, re-checks, and monitoring hospitalized patients. An alternative panel to the Prep II rotor for pre-anesthetic testing.

Analytes Include: ALT, BUN, CL-, CRE, GLU, K+, NA+, tCO2

Used For: Hospitalized patient monitoring, chronic monitoring and re-checks (for hyperadrenocorticism, diabetes, renal disease), pre-anesthetic evaluation.

Prep Profile II

Ideal For: Determining of basic health values for younger healthy patients, a basic re-check panel for some disease states and as a low cost pre-anesthetic evaluation.

Analytes Include: ALP, ALT, BUN, CRE, GLU, TP

Used For: Routine wellness screening for younger patients, lower cost pre-anesthetic screen for younger, apparently healthy patients, monitor and recheck for disease states such as diabetes mellitus and renal disease.

T4/Cholesterol

Ideal For: Routine screening of hypothyroidism in dogs and diagnostic for hyperthyroidism in cats. Titrating and monitoring patients on thyroid hormone replacement therapy or patients being treated for hyperthyroid disease.

Analytes Include: CHOL, T4

Used For: Screening for hypothyroidism in dogs, diagnosis of hyperthyroidism in cats, monitoring drug, I131 or thyroidectomy therapy.

Mammalian Liver Profile

Ideal For: Obtaining baseline values prior to administration of and concurrent with the use of NSAIDS as well as monitoring hepatic function and diagnosing liver disease.

Analytes Include: ALB, ALP, ALT, BA, BUN, CHOL, GGT, TBIL

Used For: Diagnosis and monitoring of liver disease, obtaining baseline values prior to administration of NSAIDS, and monitoring patients concurrent with the use of NSAIDS.

Avian/Reptilian Profile Plus

Ideal For: Measuring analytes that represent the most important area of concern in avian and reptilian patients.

Analytes Include: ALB, AST, BA, CA, CK, GLOB, GLU, K+, NA+, PHOS, TP, UA

Used For: Ill patient diagnostics, wellness testing, evaluate liver integrity and function, evaluate renal status, measure electrolyte status.

Equine Profile Plus

Ideal For: Routine equine checkups, wellness testing, ill patient diagnostics and pre-purchase examinations.

Analytes Include: ALB, AST, BUN, CA, CK, CRE, GLU, GGT, GLOB, K+, NA+, TBIL, tCO2, TP

Used For: Ill patient diagnostics, wellness examinations, patient monitoring, fluid therapy, and re-check examinations. Ideal for both equine ambulatory practitioners and critical care units.

Large Animal Profile

Ideal For: Health assessment, prognostic indicator and a diagnostic tool for beef and dairy cattle.

Analytes Include: ALB, ALP, AST, BUN, CA, CK, GGT, GLOB, MG, PHOS, TP

Used For: Accurate diagnosis and prognosis, aid in the choice of therapeutics.

VetCom

3240 Whipple Road
Union City, CA 94587
Telephone: 510-675-6500

Fax: 510-441-6150
vetscansales@abaxis.com
www.abaxis.com

VetCom is now available online, bringing you the latest in diagnostic insights for the entire clinic and enlightening you with a monthly case study full of interesting information and much, much more. For a printer-friendly version of VetCom to share with your colleagues, visit our website at www.abaxis.com. To sign up for additional Abaxis information please email Valerie Goodwin at valeriegoodwin@abaxis.com.

You can unsubscribe from our contact database by replying to the e-mail with the word "unsubscribe" in the subject line or simply call (800) 822-2947. For more on our Privacy Policy please visit our website at www.abaxis.com.

VetScan[®] products by
ABAXIS