




## Fast. Sensitive. Easy to Use.

The VetScan Avian Influenza Rapid Test Kit is a USDA and NPIP approved test for the detection of Type A influenza virus in chickens, turkeys, ducks and geese. This robust kit is capable of identifying influenza virus from tracheal, oropharyngeal or cloacal swabs without the need for special instruments or conditions. Within 15 minutes, this colorimetric test provides accurate results allowing for immediate decisions regarding your flock.

- Useful screening, monitoring tool
  - USDA approved
  - NPIP approved
- Highly sensitive yet highly specific
  - Component of comprehensive testing and surveillance program
  - Allows for reliable in-house screening program
  - Allows for price competitive protocol
- Speed and ease of use
  - Results easily interpreted
  - Results within 15 minutes
  - No special equipment or conditions required
  - Accuracy when used by personnel with minimal training

The Abaxis Animal Health portfolio consists of the VetScan VS2 chemistry, electrolyte and immunoassay analyzer, the VetScan VSpro, a coagulation and specialty analyzer with a combination PT/aPTT test cartridge, the VetScan i-STAT® 1, a handheld clinical analyzer, the VetScan HM2 (3-part differential) and the VetScan HM5 (5-part differential) hematology analyzer.



[www.abaxis.com](http://www.abaxis.com)

**Abaxis, Inc. - Animal Health  
Worldwide Headquarters**

3240 Whipple Road  
Union City, CA 94587  
Tel 800 822 2947  
Fax 510 441 6150