



The Piccolo Xpress delivers the vital chemistry information you need right now! **No waiting!**

Managing the health of your oncology patients is a delicate balance of treatment and monitoring. With the Piccolo Xpress in your facility, you can run comprehensive diagnostics on your patients on-site in minutes to ensure that treatments are based on the most **timely and accurate** information, which can lead to more consistent progress and a reduction in complications.

With our **Basic Metabolic Panel Plus**, you can quickly check your patient's potassium, sodium, calcium, chloride, kidney function, blood sugar, total carbon dioxide as well as magnesium and lactate dehydrogenase levels with a single panel. This information is vital to the well being of your patients and will help you treat them faster and better, so they can spend more time at home... not in your office.

Exceptionally Fast

From sample to results in 3 simple steps and approximately 12 minutes



Exceptionally Easy

Fully automated system requires no special operating skills

- Self-contained BMP Plus panel during treatment and CMP panel for follow-up patients
- Requires a very small sample - only 0.1cc whole blood, serum or plasma
- Barcoded reagent disc self calibrates with every test run
- Built-in iQC[™] eliminates daily control testing and meets CLIA and NCCLS guidelines
- Results are printed on a self-adhesive hard copy or via electronic transfer, allowing immediate access to lab-accurate chemistry results

Exceptionally Accurate

Advanced technology ensures reliable results

- Gold-standard "wet chemistry" results are measured on 9 different wavelengths
- Accuracy and precision comparable to larger, more expensive laboratory analyzers
- iQC[™] system performs continuous quality control checks on every run