

i-STAT® Level 1 Control

**LOT** 101098

 **Exp.** 2019-02-28

**CLEW: A35**

<b>CG4+ &amp; G3+</b>		<b><math>\bar{x}</math> (Mean)</b>	<b>R (Range)</b>
<b>LOT</b>	<b>D</b>		
pH		7.036	6.986 - 7.086
<b>PCO2</b>	mmHg	56.8	49.3 - 64.3
	kPa	7.57	6.57 - 8.57
<b>PO2</b>	mmHg	87	72 - 102
	kPa	11.6	9.6 - 13.6
Lactate/Lac	mmol/L	7.99	6.91 - 9.07
	mg/dL	72.0	62.3 - 81.7
	g/L	0.720	0.623 - 0.817

**CONTROL**

**IVD**



Manufactured for:  
Abbott Point of Care Inc.  
100 and 200 Abbott Park Rd.  
Abbott Park, IL  
60064 USA

**EC REP**

EMERGO EUROPE  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



748726-01

20-Dec-2017

i-STAT® Level 1 Control

**LOT** 101098

 **Exp.** 2019-02-28

**CLEW: A35**

<b>CG8+, EG7+, EG6+, G3+, Crea &amp; CG4+</b> <b>LOT A, K, L, M, N, W, Y</b>		<b><math>\bar{x}</math> (Mean)</b>	<b>R (Range)</b>
Na	mmol/L, mEq/L	122	118 - 127
K	mmol/L, mEq/L	2.8	2.5 - 3.1
iCa	mmol/L	1.48	1.36 - 1.59
	mg/dL	5.9	5.4 - 6.4
	mEq/L	3.0	2.7 - 3.2
pH		7.039	6.989 - 7.089
PCO2	mmHg	58.0	50.5 - 65.5
	kPa	7.73	6.73 - 8.73
PO2	mmHg	89	74 - 104
	kPa	11.9	9.9 - 13.9
Glucose/Glu	mg/dL	40	31 - 49
	g/L	0.40	0.31 - 0.49
	mmol/L	2.2	1.7 - 2.7
Creatinine/Crea	mg/dL	4.2	3.4 - 5.0
	umol/L	371	301 - 442
Lactate/Lac	mmol/L	8.02	6.93 - 9.10
	mg/dL	72.3	62.4 - 82.0
	g/L	0.723	0.624 - 0.820

CONTROL

IVD



Manufactured for:  
Abbott Point of Care Inc.  
100 and 200 Abbott Park Rd.  
Abbott Park, IL  
60064 USA

**EC REP**

EMERGO EUROPE  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



748726-01

20-Dec-2017

i-STAT® Level 1 Control

**LOT** 101098

 **Exp.** 2019-02-28

**CLEW: A35**

<b>CHEM8+</b> <b>LOT</b> H		<b><math>\bar{x}</math> (Mean)</b>	<b>R (Range)</b>
Na	mmol/L, mEq/L	122	117 - 126
K	mmol/L, mEq/L	2.9	2.6 - 3.2
Cl	mmol/L, mEq/L	77	72 - 82
iCa	mmol/L	1.49	1.37 - 1.60
	mg/dL	6.0	5.5 - 6.4
	mEq/L	3.0	2.7 - 3.2
TCO2	mmol/L, mEq/L	19	10 - 28
Glucose/Glu	mg/dL	41	32 - 50
	g/L	0.41	0.32 - 0.50
	mmol/L	2.3	1.8 - 2.8
Creatinine/Crea	mg/dL	4.3	3.5 - 5.1
	umol/L	380	309 - 451
BUN	mg/dL	60	52 - 68
Urea	mmol/L	21.4	18.6 - 24.3
	mg/dL	128	111 - 146
	g/L	1.28	1.11 - 1.46

**CONTROL**

**IVD**



Manufactured for:  
Abbott Point of Care Inc.  
100 and 200 Abbott Park Rd.  
Abbott Park, IL  
60064 USA

**EC REP**

EMERGO EUROPE  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



748726-01

20-Dec-2017

i-STAT® Level 1 Control

**LOT** 101098

 **Exp.** 2019-02-28

**CLEW: A35**

EC8+, 6+, EC4+, E3+ & G <b>LOT</b> F, H, J, K		$\bar{x}$ (Mean)	R (Range)
Na	mmol/L, mEq/L	122	118 - 127
K	mmol/L, mEq/L	2.8	2.5 - 3.1
Cl	mmol/L, mEq/L	75	70 - 80
pH		6.993	6.943 - 7.043
PCO2	mmHg	58.0	50.5 - 65.5
	kPa	7.73	6.73 - 8.73
Glucose/Glu	mg/dL	39	30 - 48
	g/L	0.39	0.30 - 0.48
	mmol/L	2.2	1.7 - 2.7
BUN	mg/dL	58	50 - 66
Urea	mmol/L	20.7	17.9 - 23.6
	mg/dL	124	107 - 141
	g/L	1.24	1.07 - 1.41

CONTROL

IVD



Manufactured for:  
Abbott Point of Care Inc.  
100 and 200 Abbott Park Rd.  
Abbott Park, IL  
60064 USA

**EC REP**

EMERGO EUROPE  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



748726-01

20-Dec-2017