

i-STAT® Level 3 Control

LOT 121094

 **Exp.** 2018-10-31

CLEW: A35

CG4+ & G3+		\bar{x} (Mean)	R (Range)
LOT	D		
pH		7.757	7.707 - 7.807
PCO2	mmHg	21.7	15.4 - 28.0
	kPa	2.89	2.05 - 3.73
PO2	mmHg	150	127 - 172
	kPa	20.0	16.9 - 22.9
Lactate/Lac	mmol/L	0.87	0.37 - 1.37
	mg/dL	7.8	3.3 - 12.3
	g/L	0.078	0.033 - 0.123

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748302-01

31-Oct-2017

i-STAT® Level 3 Control

LOT 121094

 **Exp.** 2018-10-31

CLEW: A35

CG8+, EG7+, EG6+, G3+, Crea & CG4+ LOT A, K, L, M, N, W, Y		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	161	156 - 166
K	mmol/L, mEq/L	6.2	5.7 - 6.7
iCa	mmol/L	0.72	0.64 - 0.80
	mg/dL	2.9	2.6 - 3.2
	mEq/L	1.4	1.3 - 1.6
pH		7.761	7.711 - 7.811
PCO2	mmHg	22.4	16.1 - 28.7
	kPa	2.99	2.15 - 3.83
PO2	mmHg	147	125 - 170
	kPa	19.6	16.7 - 22.7
Glucose/Glu	mg/dL	281	239 - 323
	g/L	2.81	2.39 - 3.23
	mmol/L	15.6	13.3 - 17.9
Creatinine/Crea	mg/dL	0.5	0.1 - 0.9
	umol/L	44	9 - 80
Lactate/Lac	mmol/L	0.87	0.37 - 1.37
	mg/dL	7.8	3.3 - 12.3
	g/L	0.078	0.033 - 0.123

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748302-01

31-Oct-2017

i-STAT® Level 3 Control

LOT 121094

 **Exp.** 2018-10-31

CLEW: A35

CHEM8+ LOT H		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	160	155 - 165
K	mmol/L, mEq/L	6.2	5.7 - 6.7
Cl	mmol/L, mEq/L	109	102 - 116
iCa	mmol/L	0.72	0.64 - 0.80
	mg/dL	2.9	2.6 - 3.2
	mEq/L	1.4	1.3 - 1.6
TCO2	mmol/L, mEq/L	33	27 - 39
Glucose/Glu	mg/dL	285	243 - 327
	g/L	2.85	2.43 - 3.27
	mmol/L	15.8	13.5 - 18.2
Creatinine/Crea	mg/dL	0.6	0.2 - 1.0
	umol/L	53	18 - 88
BUN	mg/dL	8	6 - 10
Urea	mmol/L	2.9	2.1 - 3.6
	mg/dL	17	13 - 21
	g/L	0.17	0.13 - 0.21

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748302-01

31-Oct-2017

i-STAT® Level 3 Control

LOT 121094

 **Exp.** 2018-10-31

CLEW: A35

EC8+, 6+, EC4+, E3+ & G		\bar{x} (Mean)	R (Range)
LOT	F, H, J, K		
Na	mmol/L, mEq/L	161	156 - 166
K	mmol/L, mEq/L	6.2	5.6 - 6.7
Cl	mmol/L, mEq/L	109	102 - 116
pH		7.700	7.650 - 7.750
PCO2	mmHg	22.7	16.4 - 29.0
	kPa	3.03	2.19 - 3.87
Glucose/Glu	mg/dL	279	237 - 321
	g/L	2.79	2.37 - 3.21
	mmol/L	15.5	13.2 - 17.8
BUN	mg/dL	7	5 - 10
Urea	mmol/L	2.5	1.8 - 3.6
	mg/dL	15	11 - 21
	g/L	0.15	0.11 - 0.21

CONTROL


IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748302-01

31-Oct-2017