


i-STAT® Level 3 Control

LOT 121098

 **Exp.** 2019-02-28


CLEW: A35

CG4+ & G3+		\bar{x} (Mean)	R (Range)
LOT	D		
pH		7.766	7.716 - 7.816
PCO2	mmHg	21.3	15.0 - 27.6
	kPa	2.84	2.00 - 3.68
PO2	mmHg	150	128 - 173
	kPa	20.0	17.1 - 23.1
Lactate/Lac	mmol/L	0.87	0.37 - 1.37
	mg/dL	7.8	3.3 - 12.3
	g/L	0.078	0.033 - 0.123

CONTROL


IVD



 Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748732-01

20-Dec-2017

i-STAT® Level 3 Control

LOT 121098

 **Exp.** 2019-02-28

CLEW: A35

CG8+, EG7+, EG6+, G3+, Crea & CG4+ LOT A, K, L, M, N, W, Y		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	161	156 - 167
K	mmol/L, mEq/L	6.2	5.7 - 6.8
iCa	mmol/L	0.72	0.64 - 0.80
	mg/dL	2.9	2.6 - 3.2
	mEq/L	1.4	1.3 - 1.6
pH		7.770	7.720 - 7.820
PCO2	mmHg	21.9	15.6 - 28.2
	kPa	2.92	2.08 - 3.76
PO2	mmHg	148	126 - 170
	kPa	19.7	16.8 - 22.7
Glucose/Glu	mg/dL	283	241 - 325
	g/L	2.83	2.41 - 3.25
	mmol/L	15.7	13.4 - 18.0
Creatinine/Crea	mg/dL	0.6	0.2 - 1.0
	umol/L	53	18 - 88
Lactate/Lac	mmol/L	0.85	0.35 - 1.35
	mg/dL	7.7	3.2 - 12.2
	g/L	0.077	0.032 - 0.122

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748732-01

20-Dec-2017

i-STAT® Level 3 Control

LOT 121098

 **Exp.** 2019-02-28

CLEW: A35

CHEM8+ LOT H		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	160	155 - 165
K	mmol/L, mEq/L	6.2	5.7 - 6.8
Cl	mmol/L, mEq/L	109	102 - 116
iCa	mmol/L	0.73	0.65 - 0.81
	mg/dL	2.9	2.6 - 3.2
	mEq/L	1.5	1.3 - 1.6
TCO2	mmol/L, mEq/L	33	27 - 39
Glucose/Glu	mg/dL	286	244 - 328
	g/L	2.86	2.44 - 3.28
	mmol/L	15.9	13.5 - 18.2
Creatinine/Crea	mg/dL	0.6	0.2 - 1.0
	umol/L	53	18 - 88
BUN	mg/dL	8	6 - 10
Urea	mmol/L	2.9	2.1 - 3.6
	mg/dL	17	13 - 21
	g/L	0.17	0.13 - 0.21

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748732-01

20-Dec-2017

i-STAT® Level 3 Control

LOT 121098

 **Exp.** 2019-02-28

CLEW: A35

EC8+, 6+, EC4+, E3+ & G		\bar{x} (Mean)	R (Range)
LOT	F, H, J, K		
Na	mmol/L, mEq/L	161	156 - 166
K	mmol/L, mEq/L	6.2	5.7 - 6.7
Cl	mmol/L, mEq/L	109	101 - 116
pH		7.704	7.654 - 7.754
PCO2	mmHg	21.9	15.6 - 28.2
	kPa	2.92	2.08 - 3.76
Glucose/Glu	mg/dL	280	238 - 322
	g/L	2.80	2.38 - 3.22
	mmol/L	15.5	13.2 - 17.9
BUN	mg/dL	8	5 - 10
Urea	mmol/L	2.9	1.8 - 3.6
	mg/dL	17	11 - 21
	g/L	0.17	0.11 - 0.21

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748732-01

20-Dec-2017