

i-STAT® Level 1 Control

LOT 101094

 **Exp.** 2018-10-31

CLEW: A35

CG4+ & G3+ LOT D		\bar{x} (Mean)	R (Range)
pH		7.031	6.981 - 7.081
PCO2	mmHg	56.4	48.9 - 63.9
	kPa	7.52	6.52 - 8.52
PO2	mmHg	88	73 - 103
	kPa	11.7	9.7 - 13.7
Lactate/Lac	mmol/L	7.98	6.90 - 9.06
	mg/dL	71.9	62.2 - 81.6
	g/L	0.719	0.622 - 0.816

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748298-01

31-Oct-2017

i-STAT® Level 1 Control

LOT 101094

 **Exp.** 2018-10-31

CLEW: A35

CG8+, EG7+, EG6+, G3+, Crea & CG4+ LOT A, K, L, M, N, W, Y		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	122	118 - 127
K	mmol/L, mEq/L	2.8	2.5 - 3.1
iCa	mmol/L	1.48	1.37 - 1.59
	mg/dL	5.9	5.5 - 6.4
	mEq/L	3.0	2.7 - 3.2
pH		7.035	6.985 - 7.085
PCO2	mmHg	58.3	50.8 - 65.8
	kPa	7.77	6.77 - 8.77
PO2	mmHg	90	75 - 105
	kPa	12.0	10.0 - 14.0
Glucose/Glu	mg/dL	41	32 - 50
	g/L	0.41	0.32 - 0.50
	mmol/L	2.3	1.8 - 2.8
Creatinine/Crea	mg/dL	4.2	3.4 - 5.0
	umol/L	371	301 - 442
Lactate/Lac	mmol/L	8.01	6.93 - 9.09
	mg/dL	72.2	62.4 - 81.9
	g/L	0.722	0.624 - 0.819

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748298-01

31-Oct-2017

i-STAT® Level 1 Control

LOT 101094

 **Exp.** 2018-10-31

CLEW: A35

CHEM8+ LOT H		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	121	117 - 126
K	mmol/L, mEq/L	2.8	2.5 - 3.1
Cl	mmol/L, mEq/L	77	72 - 82
iCa	mmol/L	1.49	1.38 - 1.60
	mg/dL	6.0	5.5 - 6.4
	mEq/L	3.0	2.8 - 3.2
TCO2	mmol/L, mEq/L	19	10 - 28
Glucose/Glu	mg/dL	42	33 - 51
	g/L	0.42	0.33 - 0.51
	mmol/L	2.3	1.8 - 2.8
Creatinine/Crea	mg/dL	4.2	3.4 - 5.0
	umol/L	371	301 - 442
BUN	mg/dL	60	52 - 68
Urea	mmol/L	21.4	18.6 - 24.3
	mg/dL	128	111 - 146
	g/L	1.28	1.11 - 1.46

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748298-01

31-Oct-2017

i-STAT® Level 1 Control

LOT 101094

 **Exp.** 2018-10-31

CLEW: A35

EC8+, 6+, EC4+, E3+ & G LOT F, H, J, K		\bar{x} (Mean)	R (Range)
Na	mmol/L, mEq/L	122	118 - 126
K	mmol/L, mEq/L	2.8	2.5 - 3.1
Cl	mmol/L, mEq/L	75	70 - 80
pH		6.988	6.938 - 7.038
PCO2	mmHg	58.7	51.2 - 66.2
	kPa	7.82	6.82 - 8.82
Glucose/Glu	mg/dL	40	31 - 49
	g/L	0.40	0.31 - 0.49
	mmol/L	2.2	1.7 - 2.7
BUN	mg/dL	58	50 - 66
Urea	mmol/L	20.7	17.9 - 23.6
	mg/dL	124	107 - 141
	g/L	1.24	1.07 - 1.41

CONTROL

IVD



Manufactured for:
Abbott Point of Care Inc.
100 and 200 Abbott Park Rd.
Abbott Park, IL
60064 USA

EC REP

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



748298-01

31-Oct-2017